

Product Specifications

Item #MM005



MAGID

Magid® MM005 Disposable Medical Masks

PRODUCT DESCRIPTION

Product Classification	Respiratory Protection - Disposable Medical Masks
Construction/Design	Triple-layer construction, ear loops
Composition	Non-woven material inner and outer mask, Spandex/polyester blend ear loops, galvanized iron adjustable nose clip
Color(s)	Blue/white
Unit of Measure	Each (50 each minimum order)
Packaging	50 each per box, 2000 each per case

FEATURES & BENEFITS

- Disposable medical masks shield the mouth and nose from large-particle droplets, splashes, sprays, or splatters
- Triple-layer construction
- Non-woven material
- Comfortable stretch polyester ear loops
- Adjustable metal nose clip prevents eyewear from fogging
- Please refer to General Information below for notes on proper usage

SERVICE LIFE

The service life of this product will vary depending on the application and type of care administered. This product should be replaced immediately once a critical or major defect (as defined below) is noted.

Critical Defects	Holes/tears
Major Defects	Dimensional or weight nonconformance, discoloration, other fabric defects

FUNCTIONAL PERFORMANCE

EC Production Quality Assurance System

Standard	Result
Directive 93/42/EEC on Medical Devices, Annex V (Devices in Class I in sterile conditions, sterilized systems or procedure packs)	Certified

GENERAL INFORMATION

Notes on Proper Usage: If worn properly, a medical mask is meant to help block large-particle droplets, splashes, sprays, or splatter that may contain germs (viruses and bacteria), keeping it from reaching your mouth and nose. Medical masks may also help reduce exposure of your saliva and respiratory secretions to others.

While a medical mask may be effective in blocking splashes and large-particle droplets, a face mask, by design, does not filter or block very small particles in the air that may be transmitted by coughs, sneezes, or certain medical procedures. Medical masks also do not provide complete protection from germs and other contaminants because of the loose fit between the surface of the face mask and your face.

Medical masks are not intended to be used more than once. If your mask is damaged or soiled, or if breathing through the mask becomes difficult, you should remove the mask, discard it safely, and replace it with a new one. To safely discard your mask, place it in a plastic bag and put it in the trash. Wash your hands after handling the used mask.

These masks do not supply oxygen.

Originator	Scott Daigle	Date	3/31/2020	
Revised By		Date		Revision #

Please note: Specifications are subject to change.

Neither this nor any statement made by or on behalf of Magid Glove & Safety herein is a warranty that any product is fit for a particular purpose, and Magid Glove & Safety shall not be liable to anyone with respect to this information.

© Magid Glove & Safety Mfg. Co. LLC. Magid® is a registered trademark of Magid Glove & Safety Mfg. Co. LLC. All other trademarks are the property of their respective owners. This document is the property of Magid Glove & Safety and the information contained herein is proprietary to Magid.

WD